

Quantitative determination of LDL Cholesterol in human Serum / Plasma / other body fluids. Only for *In Vitro* Diagnostic use

ORDER INFORMATION

| REF | Pack Size |
|-----------|------------|
| LDL 20 | 1 X 20 ml |
| LDL 40 | 1 X 40 ml |
| LDL 100 | 1X100 ml |
| LDL 1000 | 1X1000 ml |
| LDL 5000 | 1X5000 ml |
| LDL 10000 | 1X10000 ml |

CLINICAL SIGNIFICANCE

Lipoproteins serve to solubilise and transport cholesterol and other lipids in the bloodstream. Different lipoprotein classes have been shown to have various effects on the progression of coronary heart disease (CHD).

The measurement of LDL cholesterol (LDL-C) is a powerful predictor of CHD that provides the opportunity for early diagnosis and intervention to halt the progress of cardiovascular disease

Method

Photometric Test method. The **LDL Cholesterol** Reagent is based on a modified polyvinyl sulfonic acid (PVS) and polyethylene-glycol methyl ether (PEGME) coupled classic precipitation method with the improvements in using optimized quantities of PVS/PEGME and selected detergents. LDL, VLDL, and chylomicron (CM) react with PVS and PEGME and the reaction results in inaccessibility of LDL, VLDL and CM by cholesterol oxidase (CHOD) and cholesterol esterase (CHER), whereas HDL reacts with the enzymes. Addition of R2 containing a specific detergent releases LDL from the PVS/PEGME complex. The released LDL reacts with the enzymes to produce H₂O₂ which is quantified by the Trinder reaction.

PRINCIPLE

The LDL Direct Cholesterol assay is a homogeneous method for directly measuring LDL-C levels in serum or plasma, without the need for any off-line pretreatment or centrifugation steps.

The method is in a two reagent format and depends on the properties of a unique detergent. This detergent (Reagent 1) solubilizes only the non LDL lipoprotein particles. The cholesterol released is consumed by cholesterol esterase and cholesterol oxidase in a non color forming reaction.

A second detergent (Reagent 2) solubilizes the remaining LDL particles and a chromogenic coupler allows for color formation. The enzyme reaction with LDL-C in the presence of the coupler produces color that is proportional to the amount of LDL cholesterol present in the sample.

REAGENT

- Reagent 1 : Enzyme Reagent
- Reagent 2 : Developer Reagent
- Calibrator : (Lyophilized) Human Serum.

REAGENT PREPARATION

Reagents are ready to use as supplied. **Reconstitute the calibrator with the exact volume of deionized water as mentioned on the label.** Mix well. Allow to stand at room temperature for 30 minutes.

REAGENT STORAGE AND STABILITY

When stored at 2-8°C reagent is stable until the expiration date stated on the bottle and kit box label.

WARNING AND PRECAUTIONS

- For in vitro diagnostic use.
- Do not use components beyond the expiration date.
- Do not mix materials from different kit lot numbers.
- Exercise the normal precautions required for handling all laboratory reagents.
- The reagent contains preservative. Do not swallow. Avoid contact with skin and mucous membranes.
- For detailed information refer Material Safety Data Sheet.

WASTE MANAGEMENT

Please refer to local legal requirements.

MATERIALS REQUIRED BUT NOT PROVIDED

- NaCl solution 9 g/L
- General laboratory equipment

SAMPLE COLLECTION AND PRESERVATION

Serum or heparin plasma

It is very important to store the sample protected from light!

Stability: 1 day at 20 – 25°C

7 days at 4 – 8°C

3 months at –20°C in case of immediate freezing.

Freeze only once! Discard contaminated specimens!

ASSAY PROCEDURE

Operating Instructions

- Check reagent inventories at least daily to ensure that quantities are sufficient for the planned work load.
- Bring all reagents, Calibrator and samples to room temperature 18 - 28°C, prior to analysis.

| AUTOMATED PARAMETERS | |
|----------------------|-----------------------|
| Wavelength | 546 nm |
| Measurement | Against Reagent blank |
| Cuvette | 1 cm light path |
| Reaction Temperature | 37°C |
| Reaction Type | End Point |
| Reaction Direction | Increasing |
| Incubation | 5 Min. + 5 Min. |
| Sample Volume | 10 µl |
| Reagent I Volume | 750 µl |
| Reagent II Volume | 250 µl |
| Low Normal | 0 mg/dl |
| High Normal | 100 mg/dl |
| Linearity | 400 mg/dl |

MANUAL ASSAY PROCEDURE

Pipette into Test Tubes

| | CALIBRATOR | SAMPLE |
|--|------------|--------|
| REAGENT I | 750 µl | 750 µl |
| CALIBRATOR | 10 µl | - |
| SAMPLE | - | 10 µl |
| Mix well and incubate for 5 mins at 37°C & Immediately Add | | |
| REAGENT II | 250 µl | 250 µl |
| After 5 minutes at 37°C. Read the absorbance (Ac) for calibrator and absorbance (As) for sample. | | |

SAMPLE DILUTIONS

- This method is linear upto a concentration of 400 mg/dL.
- Dilute samples above this concentration 1:1 with 0.9% saline
- Repeat assay. Multiply the result by 2.

CALCULATION

Results are calculated, usually automatically by the instrument, as follows:

$$\text{LDL Chol (mg/dL)} = \frac{(A2-A1) \text{ of Unknown} \times \text{Calibrator value}}{(A2-A1) \text{ of Calibrator}}$$

CLIBRATORS AND CONTROLS

For the calibration of automated photometric systems the commercially available suitable multi-calibrator is recommended.

It is recommended to run a normal and a pathological control serum which is commercially available to verify the performance of the measured procedure. The value of controls should fall within the established limit. Each laboratory should establish corrective action in case of deviations in control recovery.

PERFORMANCE CHARACTERISTICS

WITHIN RUN

| Sample | Mean Concentration | SD | CV % |
|------------------------------|--------------------|------|-------|
| Randox Lipid Control Level 1 | 85.51 | 0.81 | 0.95% |
| Randox Lipid Control Level 2 | 138.00 | 0.98 | 0.71% |

RUN TO RUN

| Sample | Mean Concentration | SD | CV % |
|------------------------------|--------------------|------|-------|
| Randox Lipid Control Level 1 | 85.36 | 0.26 | 0.30% |
| Randox Lipid Control Level 2 | 139.36 | 0.45 | 0.32% |

LINEARITY

The method is linear upto a concentration of 400mg/dL. Dilute samples above this concentration 1:1 with 0.9% saline solution and repeat assay. Multiply the result by 2.

Limit of detection: The limit of detection for LDL Cholesterol is 2 mg/dL.

METHOD COMPARISON

A comparison of LDL Cholesterol with a commercially available assay (x) using 20 samples gave following results: $R^2 = 0.9900$

REFERENCE VALUES

| | |
|-----------------------|---------|
| Optimal | <100 |
| Near or above optimal | 100-129 |
| Borderline high | 130-159 |
| High | 160-189 |
| Very high | >190 |

The reference values are to be considered as indicative only. Every laboratory should establish its own normal range.

LIMITATION OF THE PROCEDURE

- For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.






INTERFERENCE

- Hemoglobin: No interference found upto 500 mg/dL.
- Bilirubin: No interference found upto Bilirubin (free) 50 mg/dL, Bilirubin (Conjugated) 40m/dL.
- Ascorbic Acid: No interference found upto 50 mg/dL.
- These characteristics have been obtained using an automatic analyzer. Results may vary if a different instrument or a manual procedure is used.

BIBLIOGRAPHY

Friedewald, W.T., Levy, R.I., and Fredrickson, D.S., Estimation of the low density lipoprotein cholesterol in plasma without use of the preparative ultracentrifuge, Clin. Chem. 18, 1972, p. 449 – 502

GLOSSARY OF SYMBOL

| | |
|---|-----------------------------|
|  | Consult Instruction for Use |
|  | Catalog Number |
|  | Store between |
|  | Manufacturer |
|  | Keep away from sunlight |



Paramcare Life Sciences Private Limited, G/F-12/13, Evershine-2,
Survey No. 307/3/1, Balitha N.H No 48, Vapi, Valsad, Gujarat, 396191.

Quanton Biolife Sciences Private Limited
Anand Mangal Apartment, Behind Axis Bank,
Dak Bunglow Road, Ghatsila, East Singhbhum
Jharkhand - 832303 India
quantoncare@qblsci.com
www.quantonbiolifesciences.com